Acorn Stairlifts is essential and open for you! Call us Today!

(/covid)

Outdoor Stairlifts: Costs & Installation

Open up and enjoy your outdoor space with Acorn Stairlifts.



(/assets/img/outdoor/zoom/acorn-outdoor-stair-lift-uk.jpg)

Do you have a steep drive that is difficult to walk up? Acorn Stairlifts do not need to be limited to your indoor stairs; our outdoor stairlifts are the perfect mobility aid to give you access to those outdoor spaces you love or have limited access to.

Our 130 outdoor stairlift comes with all the features that can be found on our **straight stairlift** (/stairlift/stairlift-straight-staircase). For extra protection, it comes with a weather resistant outdock cover.

View Technical Details

Why not request your free brochure today?

Safety and Convenience

- 1. Folding foot rest
- 2. Seat belt
- 3. Diagnostic display
- 4. Backup battery
- 5. Remote control

Outdoor Stairlift Features

- 1. Fits quickly and easily to the stairs, not the wall.
- 2. Stairlift installation requires no major structural changes.
- 3. Our stairlift stops automatically if it meets an obstruction.
- 4. Slimline and folds away conveniently when not in use.
- 5. Optional hinged rail designed to eliminate any obstruction or trip hazard.

Outdoor Stairlift Installation

All Acorn stairlifts fit directly to the staircase, not the wall, so installation is quick and mess-free.


FREE Outdoor Stairlift Quote

Our comprehensive **home survey (/contact/free-stairlifts-survey)** is completely free and you're under no obligation to buy!


FREE Stairlift Home Survey FREE Stairlift Quote NO Pressure Selling NO Need for Repairs or Redecoration

Contact us for FREE no obligation advice and speak to one of our friendly experts.

FREEPHONE: 0808 223 4871 (tel:08082234871)



PDF 🔏 ACORN OUTDOOR STAIRLIFT 130 - USER MANUAL (../Acorn-Stairlifts-User-Manuals-130.pdf)



Angles

Standard Indoor 130

Short seat post	40-55°
Standard seat post	29-55°
Tall seat post	23-55°

Power Swivel 130

Standard seat post	36-55°
Tall seat post	23-55°

Øutdoor 130

Standard seat post

27-55°

Maximum Rail Length

10	0 metres 394 inches			
Dir	Dimensions mm inch			
А	A Floor to top of footrest		95	3.75
В	Top of footrest to top of seat (Short) (Tall)		445 (410) (505)	17.5 (16.25) (19.75)
С	Top of seat to top of arms		250	9.75
D	Top of arms to top of seat back		190	17.5
Е	Width between armrests (Wide arm version)		445 (495)	17.5 (19.5)
F	Overall width		605	23.75
G	Overall height to floor		980	38.5
Н	Length of arms from seat back		470	18.5
I	Length of seat base from seat back		435	17
J	Length of footrest		325	12.75
К	Width of footrest		310	12.25
L	Front of footrest to stringer		570	22.5
Μ	Back of seat to stringer		13	0.5
Ν	Folded width		310	12.25
0	Minimum width required to swivel at top		670	26.25
Ρ	Distance from front of rail to stringer (Wide ar	m version)	215 (240)	8.5 (9.5)

This stairlift is independently tested and manufactured to the following standards:

- BS EN 81-40:2020 European safety standard 'Specification for powered stairlifts'.
- CE mark Acorn Stairlift meets the requirements of the applicable European Directive(s).
- ISO 9386:2-2000 International standard for stairlifts.

Specifications Motor output speed No greater than 0.15 metres per second

(/) Specifications			
Method of drive	Rack and Pinion		
Power supply	24v DC batteries		
Motor Power	0.25KW		
Maximum capacity	127kg		
Maximum capacity HD version	159kg		
Power supply input range	50-60Hz, 120-240V AC		

* NOTE: Acorn, in keeping with its policy of continual development, reserves the right to change specification without notice. All measurements are approximate.



Get Your FREE Quote Now

First Name

Last Name

1/4

NEXT

Why choose the world's most innovative stairlift?



Fits quickly and easily to the stairs, not the wall

Requires no structural changes to your home Affordable and reliable with low running costs

Slimline and folds away conveniently when not in use

Learn what our happy stairlift owners say...

If you want to know how an Acorn Stairlift can transform your life, don't take our word for it — see what our customers have to say...

81 year old Margaret Sergison found she could no longer tackle her staircase, often finding herself immobile. A dangerous fall down the stairs forced her eldest son, Michael, to take action.

He contacted Acorn Stairlifts who installed his mother's stairlift just one week later. "The difference is a blessing," says Mrs Sergison.

"I simply can't explain how much the Acorn Stairlift has helped me. Most importantly, I can stay in the home I love, which is my life."

READ MORE CUSTOMER STORIES (/ABOUT/ACORN-STAIRLIFTS-REVIEWS)

Acorn Information

Our Stairlifts

Follow Acorn

- About Acorn Stairlifts (/about/about-acorn-stairlifts)
- Existing Customers (/about/acorn-stairlift-existing-customer)
- Quality & Testing (/stairlift/stairlift-quality)
- Complaints Procedure (/contact/complaints)
- Dr. Hilary Jones (/about/dr-hilary-jones)
- Our Accreditations (/about/acorn-stairlifts-awards)
- Contact Us (/contact/contact-us)

- (/) Privacy Policy (/info/privacy-policy)
 - Cookie Policy (/info/cookies-policy)
 - Anti-Slavery Policy (/info/anti-slavery-policy)
 - Gender Pay Gap Reports (/info/gender-pay-gap)
 - Careers (/about/acorn-stairlifts-careers)
 - Our Locations (/info/our-locations)
 - Sitemap (/sitemap.xml)
 - Terms & Conditions (/info/TCs)
 - Safeguarding Policy (/info/safeguarding-policy)

Acorn Mobility Services Limited Telecom House Millennium Business Park Steeton West Yorkshire, UK BD20 6RB

FREEPHONE: 0808 223 4871 (tel:08082234871)

Registered in England No. 02593771

VAT No. 607057749

General Enquires:

customerrelations@acornstairlifts.co.uk

(mailto:customerrelations@acornstairlifts.co.uk)

(/)

All Rights Reserved © 2021 Acorn Mobility Services Ltd